

# Welcome to *activforce 2*.

## In this kit, you will find:

- The Activeforce 2 device
- The standard shell attachment
- The cylindrical shell for measuring a wrist, head and other angular body parts
- A hand strap for securely attaching the device to your hand
- A leg strap to secure the device to a patient's leg or other body part
- A web belt for measuring pull force
- A belt loop plate to secure the straps or belt to the device

## To Get Started:

- 1 Download the Activforce 2 app




- 2 Create an account or log into an existing account
- 3 Turn your device on by holding the blue button on the side of your device down until you see a blue light and tap the "Connect" button on the app
- 4 Complete the onboard questions
- 5 Start testing


For How-to guides and other information, visit [activbody.com/docs](http://activbody.com/docs) or email us at [support@activbody.com](mailto:support@activbody.com)

# Activforce 2 Accessories


Shell Attachment



Cylindrical Attachment




Belt Loop Plate



Hand Strap



Leg Strap



Belt

